











B.LV5.F.GL.FC.S

Biosensor array

For bioanalytical applications a flow-through sensor


A product of Jobst Technologies


Benefits & Characteristics

- Enzyme-based amperometric measurement
- Excellent long-term stability
- Stable in continuous monitoring and analyzer mode
- Reference, counter and blank electrodes on-chip
- Integrated flow cell (various volumes)
- For industrial applications

- Suitable for multiparametric measurements
- Outstanding reliability
- Fast response time
- Gamma and beta sterilization compatible
- Suitable for flow-through applications

Illustration¹⁾





1) For actual size, see dimensions

Technical Data²⁾

Dimensions - chip (L x W x H in mm):	10 x 4 x 0.4 (1.23 mm with flow cell)	
Dimensions - flow cell (V _{INTERNAL}) ³⁾ :	0.3 μL to 1 μL	~ 8 mm x 600 µm (0.15, 0.5 mm tubes)
Fluid connection (Ø _{INNER}):	0.5 mm (0.15 mm available upon request)	
Working electrode:	Platinum covered with enzyme membrane	
Blank electrode:	For background compensation	
Reference electrode:	Silver / silver chloride	
Counter electrode:	Platinum	
Measurable analytes ³⁾ :	Glucose, Lactate, Glutamine ⁴⁾ , Glutamate ⁴⁾	
Operating measurement range at +25 °C (for extended measurement ranges please contact IST AG) ⁵⁾ :	Glucose:	0.05 mM to 25 mM / 0.01 to 4.5 g/l
	Lactate:	0.02 mM to 15 mM / 0.002 to 1.5 g/l
	Glutamine:	0.04 mM to 20 mM / 0.006 to 3g/l
	Glutamate:	0.04 mM to 10 mM / 0.006 to 1.5 g/l













Sterilization:	Beta or gamma sterilization is possible (25 kGy typical)
	Sensitivity increases upon radiation; lactate sensors life
	time is reduced
	Initial bioburden < 1 cfu per sensor

	initial biobarden virieta per sensor
Time to first measurement after storage at +37 °C:	~ 15 min
Time response (t_{90} %) at +37 °C ⁵⁾ :	< 25 s
Temperature influence:	~ 3.5 %/°C
Storage conditions:	+4 °C to +35 °C, desiccated
Shelf life:	> 6 months (at recommended storage conditions)
Operational life time ⁶⁾ :	> 2 weeks continuous operation> 4 weeks in analyzer mode

²⁾ All performance data was obtained using an acetate buffer. For more information about medium composition please contact IST AG

Flow cell pressure drop:

Product Photo



Sensor only



 \sim 700 μ l/min/bar (with a 1 μ l flow cell)

Packaged sensor

³⁾ Larger cell volumes available upon request

⁴⁾ Upon request

⁵⁾ Parameters may vary depending on the medium type

⁶⁾ Limited by lactate sensor; may decrease by higher analyte concentration and / or less dissolved oxygen in medium














Pin Assignment⁷⁾




1	2	3	4
Counter	Reference	Lactate 2	Glucose 2
5	6	7	8
Blank 2	Lactate 1	Glucose 1	Blank 1

⁷⁾ Customer specific pin assignment upon request

Packaged Sensor

Innovative Sensor Technology IST AG offers a highly robust packaged version of the LV5 sensor for easy mounting and testing.



Electrical connection:	Card connector compatible with MEC6-RA connectors
	(40 pins in two rows, pitch: 0.635 mm)

Fluidic connection: 1 x female luer lock, 1 x male luer lock













Order Information - Sensor only

Description: Item number: Former main reference: B.LV5.F.GL.FC 105117 390.00129

Order Information - Packaged Sensor

Description: Item number: Former main reference: B.LV5.F.GL.FC.S 105128 390.00154

Disclaimer

Not for medical, diagnostics and use on humans. For evaluation use only. For more information contact IST AG.



